

RAB40B Antibody



PACO11662

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RAB40B

Uniprot

Q12829

Synonyms:

RAB40B, member RAS oncogene family; RAB40B; FLJ42385; RAR; SEC4L ;

Immunogen:

Human RAB40B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A