

RABEP1 Antibody



PACO11676

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RABEP1

Uniprot

Q15276

Synonyms:

rabaptin, RAB GTPase binding effector protein 1; RABEP1; RAB5EP; RABPT5 ;

Immunogen:

Human RABEP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A