

RABEP2 Antibody



PACO11677

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

RABEP2

Uniprot

Q9H5N1

Synonyms:

rabaptin, RAB GTPase binding effector protein 2; RABEP2; FLJ23282; FRA ;

Immunogen:

Human RABEP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A