

# RABGAP1L Antibody



PACO11679

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RABGAP1L

**Uniprot**

B7ZAP0

**Synonyms:**

RAB GTPase activating protein 1-like; RABGAP1L; DKFZp686E1450; HHL; TBC1D18 ;

**Immunogen:**

Human RABGAP1L.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A