

# RABGEF1 Antibody



PACO11680

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RABGEF1

**Uniprot**

Q9UJ41

**Synonyms:**

RAB guanine nucleotide exchange factor (GEF) 1; RABGEF1; FLJ32302; RABEX5; rabex-5 ;

**Immunogen:**

Human RABGEF1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A