

# RABL2A Antibody



PACO11682

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RABL2A

**Uniprot**

Q9UBK7

**Synonyms:**

RAB, member of RAS oncogene family-like 2A; RABL2A; FLJ78724; MGC117180 ;

**Immunogen:**

Human RABL2A.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A