

# RAD51B Antibody



PACO11696

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RAD51B

**Uniprot**

O15315

**Synonyms:**

RAD51-like 1 (*S. cerevisiae*); RAD51L1; MGC34245; R51H2; RAD51B; REC2; hREC2 ;

**Immunogen:**

Human RAD51L1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A