

RAI1 Antibody



PACO11699

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RAI1

Uniprot

Q7Z5J4

Synonyms:

retinoic acid, induced 1; RAI1; DKFZp434A139; KIAA1820; MGC12824; SMCR; SMS ;

Immunogen:

Human RAI1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A