

# RANBP3L Antibody



PACO11710

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RANBP3L

**Uniprot**

Q86VV4

**Synonyms:**

RAN binding protein 3-like; RANBP3L; FLJ25422 ;

**Immunogen:**

Human RANBP3L

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A