

RANGRF Antibody



PACO11714

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RANGRF

Uniprot

Q9HD47

Synonyms:

RAN guanine nucleotide release factor; RANGRF; DKFZp686F02139; HSPC165; HSPC236; MGC110973; MOG1; RANGNRF ;

Immunogen:

Human RANGRF.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A