

RASA3 Antibody



PACO11729

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RASA3

Uniprot

Q14644

Synonyms:

RAS p21 protein activator 3; RASA3; GAP1IP4BP; GAPIII; MGC46517; MGC47588 ;

Immunogen:

Human RASA3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A