

# RASGEF1A Antibody



PACO11731

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RASGEF1A

**Uniprot**

Q8N9B8

**Synonyms:**

RasGEF domain family, member 1A; RASGEF1A; CG4853; FLJ37817 ;

**Immunogen:**

Human RASGEF1A.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A