

# RASGEF1B Antibody



PACO11732

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

RASGEF1B

**Uniprot**

Q0VAM2

**Synonyms:**

RasGEF domain family, member 1B; RASGEF1B; FLJ31695; GPIG4; MGC46251 ;

**Immunogen:**

Human RASGEF1B.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A