

RASSF7 Antibody



PACO11737

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RASSF7

Uniprot

Q02833

Synonyms:

Ras association (RalGDS/AF-6) domain family (N-terminal) member 7; RASSF7; C11orf13; HRAS1; HRC1; MGC126069; MGC126070 ;

Immunogen:

Human RASSF7.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A