

RASSF9 Antibody



PACO11739

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RASSF9

Uniprot

O75901

Synonyms:

Ras association (RalGDS/AF-6) domain family (N-terminal) member 9; RASSF9; P-CIP1; PAMCI ;

Immunogen:

Human RASSF9.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A