

RBBP4 Antibody



PACO11742

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RBBP4

Uniprot

Q09028

Synonyms:

retinoblastoma binding protein 4; RBBP4; NURF55; RBAP48 ; CAF-I p48

Immunogen:

Human RBBP4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A