

RBM39 Antibody



PACO11759

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RBM39

Uniprot

Q14498

Synonyms:

RNA binding motif protein 39; RBM39; CAPER; CAPERalpha; CC1.3; DKFZp781C0423; FLJ44170; HCC1; RNPC2; fSAP59 ;

Immunogen:

Human RBM39.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A