

RBM41 Antibody



PACO11761

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

RBM41

Uniprot

Q96IZ5

Synonyms:

RNA binding motif protein 41; RBM41; FLJ11016; FLJ13670; bB383K5.1 ;

Immunogen:

Human RBM41.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A