

RBM42 Antibody



PACO11762

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RBM42

Uniprot

Q9BTD8

Synonyms:

RNA binding motif protein 42; RBM42; MGC10433 ;

Immunogen:

Human RBM42.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A