

RBMS2 Antibody



PACO11767

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RBMS2

Uniprot

Q15434

Synonyms:

RNA binding motif, single stranded interacting protein 2; RBMS2; FLJ39093; FLJ40023; FLJ43262; SCR3 ;

Immunogen:

Human RBMS2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A