

RCAN2 Antibody



PACO11777

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RCAN2

Uniprot

Q14206

Synonyms:

regulator of calcineurin 2; RCAN2; CSP2; DSCR1L1; MCIP2; RCN2; ZAKI-4; ZAKI4; hRCN2 ;

Immunogen:

Human RCAN2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A