

RDH11 Antibody



PACO11785

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RDH11

Uniprot

Q8TC12

Synonyms:

retinol dehydrogenase 11 (all-trans/9-cis/11-cis); RDH11; ARSDR1; CGI-82; FLJ32633; HCBP12; MDT1; PSDR1; RALR1; SCALD; SDR7C1 ;

Immunogen:

Human RDH11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A