

RDH13 Antibody



PACO11787

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

RDH13

Uniprot

Q8NBN7

Synonyms:

retinol dehydrogenase 13 (all-trans/9-cis); RDH13; DKFZp313H0740; FLJ39479; SDR7C3 ;

Immunogen:

Human RDH13.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A