

RDH14 Antibody



PACO11788

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

RDH14

Uniprot

Q9HBH5

Synonyms:

retinol dehydrogenase 14 (all-trans/9-cis/11-cis); RDH14; PAN2; SDR7C4 ;

Immunogen:

Human RDH14.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A