

RHPN1 Antibody



PACO11838

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RHPN1

Uniprot

Q8TCX5

Synonyms:

rhophilin, Rho GTPase binding protein 1; RHPN1; KIAA1929; ODF5; RHOPHILIN; RHPN ;

Immunogen:

Human RHPN1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A