

RIBC2 Antibody



PACO11839

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

RIBC2

Uniprot

Q9H4K1

Synonyms:

RIB43A domain with coiled-coils 2; RIBC2; C22orf11 ;

Immunogen:

Human RIBC2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A