

# RIC3 Antibody



PACO11840

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RIC3

**Uniprot**

Q7Z5B4

**Synonyms:**

resistance to inhibitors of cholinesterase 3 homolog (C. elegans); RIC3; AYST720; FLJ11608; PRO1385 ;

**Immunogen:**

Human RIC3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A