

RIC8B Antibody



PACO11842

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

RIC8B

Uniprot

Q9NVN3

Synonyms:

resistance to inhibitors of cholinesterase 8 homolog B (C. elegans); RIC8B; FLJ10620; MGC39476; RIC8; hSyn ;

Immunogen:

Human RIC8B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A