

RIMS3 Antibody



PACO11843

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RIMS3

Uniprot

Q9UJD0

Synonyms:

regulating synaptic membrane exocytosis 3; RIMS3; KIAA0237; NIM3; RIM3 ;

Immunogen:

Human RIMS3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A