

RNASE13 Antibody



PACO11857

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RNASE13

Uniprot

Q5GAN3

Synonyms:

ribonuclease, RNase A family, 13 (non-active);

Immunogen:

Human RNASE13.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A