

# RNF212 Antibody



PACO11873

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RNF212

**Uniprot**

Q495C1

**Synonyms:**

ring finger protein 212; RNF212; FLJ38841; MGC120227; MGC120228; ZHP3 ;

**Immunogen:**

Human RNF212.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A