

PTBP3 Antibody



PACO11886

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

PTBP3

Uniprot

O95758

Synonyms:

ROD1 regulator of differentiation 1 (S. pombe); ROD1; DKFZp781I1117; PTBP3 ;

Immunogen:

Human ROD1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A