

RPAP1 Antibody



PACO11897

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RPAP1

Uniprot

Q9BWH6

Synonyms:

RNA polymerase II associated protein 1; RPAP1; DKFZp727M111; FLJ12732; KIAA1403; MGC858 ;

Immunogen:

Human RPAP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A