

RPL32 Antibody



PACO11925

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RPL32

Uniprot

P62910

Synonyms:

ribosomal protein L32; RPL32; PP9932 ;

Immunogen:

Human RPL32.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A