

RPL7L1 Antibody



PACO11933

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

RPL7L1

Uniprot

Q6DKI1

Synonyms:

ribosomal protein L7-like 1; RPL7L1; MGC62004; dJ475N16.4 ;

Immunogen:

Human RPL7L1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A