

RPP40 Antibody



PACO11940

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

RPP40

Uniprot

O75818

Synonyms:

ribonuclease P/MRP 40kDa subunit; RPP40; RNASEP1; bA428J1.3 ;

Immunogen:

Human RPP40.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A