

# RPS15 Antibody



PACO11945

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

RPS15

**Uniprot**

P62841

**Synonyms:**

ribosomal protein S15; RPS15; MGC111130; RIG ;

**Immunogen:**

Human RPS15.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A