

RRAGB Antibody



PACO11965

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RRAGB

Uniprot

Q5VZM2

Synonyms:

Ras-related GTP binding B; RRAGB; RAGB; bA465E19.1 ;

Immunogen:

Human RRAGB.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A