

RSAD2 Antibody



PACO11973

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

Protein Background:

Gene ID:

RSAD2

Uniprot

Q8WXG1

Synonyms:

361aa,42kd;

Immunogen:

Human RSAD2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20oC, Avoid freeze / thaw cycles.

Product Images

N/A

N/A