

RTN1 Antibody



PACO11979

Size:

100ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Form:

Liquid

Product Type:

Polyclonal Antibody

Immunogen Species:

Human

Purification Method:

Antigen Affinity purified

Gene ID:

RTN1

Uniprot

Q16799

Synonyms:

reticulon 1; RTN1; MGC133250; NSP ;

Immunogen:

Human RTN1.

Storage Buffer:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C. Avoid freeze / thawcycles.

Storage:

Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

