

RUNDC3A Antibody



PACO11987

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

RUNDC3A

Uniprot

Q59EK9

Synonyms:

RUN domain containing 3A; RUNDC3A; RAP2IP; RPIP8 ;

Immunogen:

Human RUNDC3A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A