

RUNDC3B Antibody



PACO11988

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

RUNDC3B

Uniprot

Q96NL0

Synonyms:

RUN domain containing 3B; RUNDC3B; FLJ30671; MGC26655; RPIB9; RPIP9 ;

Immunogen:

Human RUNDC3B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A