

S100A14 Antibody



PACO12003

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

S100A14

Uniprot

Q9HCY8

Synonyms:

S100 calcium binding protein A14; S100A14; BCMP84; S100A15 ;

Immunogen:

Human S100A14.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A