

SAMD4A Antibody



PACO12019

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SAMD4A

Uniprot

Q9UPU9

Synonyms:

sterile alpha motif domain containing 4A; SAMD4A; DKFZp434H0350; KIAA1053; SAMD4; SMG; SMGA; Smaug; Smaug1 ;

Immunogen:

Human SAMD4A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A