

SAMD4B Antibody



PACO12020

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

SAMD4B

Uniprot

Q5PRF9

Synonyms:

sterile alpha motif domain containing 4B; SAMD4B; FLJ10211; MGC99832; SMGB ;

Immunogen:

Human SAMD4B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A