

SAR1A Antibody



PACO12027

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

SAR1A

Uniprot

Q9NR31

Synonyms:

SAR1 homolog A (*S. cerevisiae*); SAR1A; SAR1; SARA1; Sara; masra2 ;

Immunogen:

Human SAR1A.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A