

SCEL Antibody



PACO12043

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SCEL

Uniprot

O95171

Synonyms:

sciellin; SCEL; FLJ21667; MGC22531 ;

Immunogen:

Human SCEL.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A