

# SCLT1 Antibody



PACO12050

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SCLT1

**Uniprot**

Q96NL6

**Synonyms:**

sodium channel and clathrin linker 1; SCLT1; CAP1A; FLJ30655; FLJ36042; MGC70542; hCAP-1A ;

**Immunogen:**

Human SCLT1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A