

SCNM1 Antibody



PACO12056

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

SCNM1

Uniprot

Q9BWG6

Synonyms:

sodium channel modifier 1; SCNM1; MGC3180 ;

Immunogen:

Human SCNM1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A