

SCYL3 Antibody



PACO12065

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

SCYL3

Uniprot

Q8IZE3

Synonyms:

SCY1-like 3 (S. cerevisiae); SCYL3; PACE-1; PACE1; RP1-97P20.2 ;

Immunogen:

Human SCYL3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A